

MW43

NEW

# "AbdoAbby" Abdominal Examination Simulator

**True-to-life tactile representation of abdominal organs to cultivate the sense of touch and serve as a benchmark during palpation.**

## Product Supervision

Takahiro Amano, M.D.,  
Senior Advisor for Academic Affairs and Professor  
International University of Health and Welfare,  
Honorary Director Sannno Medical Center, Japan


Kyoko Arai  
RN, PHN, PhD, Professor,  
School of Nursing, University of Shizuoka, Japan



Interchangeable organs



Supple, resilient and delicate texture



Respiratory motion

## FEATURES

- 1 | Supple, resilient and delicate textures which allow both shallow and deep palpation
- 2 | Interchangeable organs
- 3 | Respiratory motion
- 4 | Abdominal auscultation

## SKILLS

- | Visual Inspection
- | Auscultation
  - renal artery
  - abdominal aorta
  - iliac artery
  - bowel sounds
- | Palpation
  - liver, spleen, kidney, uterus and lymphoma
- | Percussion
  - liver, spleen and rib

## KEY FEATURES

### Bony landmarks


- | Pelvis
- | Ribs
- | Lower rib cage
- | Costal margin
- | Xiphisternum
- | Pubic crest
- | Anterior superior iliac spines

### Liver

- | Normal
- | Chronic hepatitis
- | Early late cirrhosis
- | Cirrhosis

### Inguinal area

- | Normal
- | Lymphoma



### Spleen

- | Normal
- | Enlarged
- | Slightly enlarged

### Kidney

- | Normal
- | Hydronephrosis
- | Polycystic kidney

### Uterus

- | Normal
- | Myoma

### Bruits and sounds

- Renal artery
- | Abdominal aorta
- | Iliac artery
- | Bowel sounds

## DESCRIPTIONS

### SPECIFICATIONS

Size: 40x80x20cm / 15.7x31.5x7.9 inch

Weight: 14kg / 30.9lbs

### MATERIALS

Soft resin  
Latex free

### SET INCLUDES

- |   |   |
|---|---|
| 1 torso unit  | 2 inguinal lymph nodes (normal, lymphoma) |
| 1 talcum powder   | 2 uterus (normal, fibroid)                |
| 1 set of the respiratory motion control unit                  | 1 simulated stethoscope                   |
| 2 kidneys (hydronephrosis, cysts)                             | 1 controller                              |
| 2 spleens (slightly enlarged, enlarged)                       | 1 AC adapter                              |
| 4 livers (precirrhosis, cirrhosis, chronic hepatitis, normal) | 3 storage cases                           |
|   | 1 instruction manual                      |