



Nasopharyngeal Swab Collection Simulator

For training to healthcare professionals who conduct virus testing for the first time





FEATURES

- 1 | When swab reaches the right position in the nasal cavity, the swab will be pigmented in blue
- 2 | Ways to use tongue depressor with swabs can be practiced.
- 3 | The shape of nasal cavity between right and left are different.

SKILLS

Training to healthcare professionals who conduct virus testing for the first time.



DESCRIPTIONS

SPECIFICATIONS

Size: W42xD23xH38cm/W16.8x D9.2x H15.2in Weight: 2.15 kg / 4.73 lbs

MATERIALS

Soft resin, Hard resin
Latex free

SET INCLUDES

- 1 main model
- 1 cross sectional model
- 1 syringe (30ml)
- 1 disposable cup
- 10 swabs for coloring
- 1 manual
- 1 storage case

