Strap-on Catheterization / Enema Simulator

MW54 | 11435-000 (female) MW53 | 11434-000 (male) MW55 | 11436-000 (value set)



Skills

- Cleansing and aseptic technique
- Catheter insertion and placement
- Confirmation of urine discharge
- Catheter removal
- Enema

Reality in skills training combined with scenarios for communication.

Features

- Soft female genitalia permits holding the labia in open position throughout the procedure from sterilization to the placement of the catheter.
- Anatomically correct structure facilitates acquiring body sense along with the procedure.
- Compatible with balloon catheters. Water (simulated urine) flows out when appropriate length of a catheter is inserted.

Specifications

Life-sized genitalia model

W8x D21 x H18cm, 0.5kg (male) W8x D21 x H9cm, 0.4kg (female)

Set Includes

1 genitalia model 1 reservoir bag 1 drainage tube 1 strap band 1 shorts 1 lubricant

1 catheter for maintenance 1 instruction manual



The catheter can be inserted when the penis is held at the right angle.



Realistic resistance to the catheter



Injection of enema fluid

■True-to-life softness

Labia can be kept open through the catheterization procedures with the nondominal hand.





lubricant

