Instruction Manual

PH-3 Angiographic CT Head Phantom ACS

Simulated arteries for CT scanning with contrast medium (41309-100) Simulated arteries for Angiography scanning with contrast medium (41309-200) Simulated arteries for Multi energy CT (41903-300)

Specifications

a. CT head phantom · · · · · · · · · 1
Life-size adult phantom weight: 6.5kg height: approx. 33cm material: Urethane base resin, Epoxy base resin
Sample X-ray images · · · · · · · · · · · · 1
Carrying case · · · · · · · · · · · · · · · · · · ·

Instruction manual1







Introduction

- This ACS head phantom is designed for image evaluation in CT angiography scanning and can be used as educational tool for medical image interpretation.
- Three kinds of head phantoms with different density of arteries are available : CT type
- (41309-100) / Angiography type (41309-200)/MECT type (41903-300)

• Features

- The head phantom includes a synthetic skull, cervical vertebrae (C1-C7), brain and arteries, along with arteries with contrast medium embedded in the left half of the head.
- Simulated arteries cover left anterior cerebral arteries and middle cerebral artery.
- Diameters of simulated arteries range from approx.0.5 mm to 4.0 mm.
- State-of-the-art synthetic skull. (Substantia spongiosa is not reproduced.)

▲ CAUTION	
 Handle with care The phantom consists of special compositions of resin. Please handle them with the utmost care at all times. 	• Even if color on its surface might change over time, this does not affect the quality of its performance.
 Storage Store model at room temperature, away from shock, heat, moisture and direct sunlight. 	 Don't mark on the models with pen or leave any printed materials in contact with their surface. Ink marks on the models are not removable.

For inquiries and service, please contact your distributor or KYOTO KAGAKU CO.,LTD.
 KYOTO KAGAKU CONTACT & ORDERING INFORMATION


KYOTO KAGAKU co., ltd

URL http://www.kyotokagaku.com e-mail erw-kyoto@kyotokagaku.co.jp

Worldwide Inquiries & Ordering **Kyotokagaku Head Office and Factories:** 15 Kitanekoya-cho, Fushimi-ku, Kyoto, 612-8388, JAPAN Tel: +81-75-605-2510 Fax: +81-75-605-2519

All American regions **KYOTOKAGAKU AMERICA INC.** 3109 Lomita Boulevard, Torrance, CA 90505-5108, USA Tel: 1-310-325-8860 Fax: 1-310-325-8867

2021.08



%3 is the value of calculation value.