41925-000 PH-60

# Tough Whole Body Phantom "PBU-90 RUGGED"





New phantom material is designed for rough handling, improved durability and less maintenance











## **FEATURES**

| Radiology absorption and HU number approximate to human bodv.

| Main joints have close-to human articulation

| Phantom can be disassembled into 10 individual parts

## **APPLICATIONS**

Plain X-ray Basic patient positioning I Basic CT

PBU-90 allows training scenarios that involve rough handling. Details of the bones in the hands and the feet are simplified compared to those of PBU-50/60 for increased durability.

#### **ANATOMY**

skull / spine / clavicles / scapulae / ribs / sternum / coxal bones / lungs with pulmonary vessels / trachea (up to primary bronchi) / heart / liver with portal and hepatic veins / kidneys / humerus / antebrachial bone / bones of hand / femur / patella / lower leg bone / bones of foot

## **DESCRIPTIONS**

1 flat head screwdriver

#### SET INCLUDES

1 whole body phantom 1 set of sample X-ray data (DVD) head supporter manual

**SPECIFICATIONS** 

Packing size: Phantom height: W92 x D57 x H38 cm/W36 x D22 x H15 in W90 x D63 x H22 cm/W35 x D25 x H8.7 in 165 cm W89 x D57 x H16 cm/W35 x D22 x H6.3 in 65 in

Phantom weight: Packing weight: 50 kg / 110 lb 80 kg / 176 lb

## MATERIALS

Soft tissue: urethane based resin (density: 1.12) \*Phantom has no metal parts or liquid structure

#### OPTIONAL PARTS

41363-070 storage cases (consist of 2 boxes) body plates for PH-2/2B (BMI 32) 41350-200-16 body plates for PH-2/2B (BMI 40) 41350-200-17

Storage case P.12 ▶



