US-10 US-10a US-10b




Practice both transvaginal and transabdominal ultrasound procedures in one

New unit released includes normal, pathological and ectopic pregnancies.

Product Supervision


Charlotte Henningsen, MS, RT(R), RDMS, RVT, FSDMS, FAIUM Chair & Professor - Sonography Department Adventist University of Health Sciences



- 1 Compatible with actual ultrasound probe
- 2 Three types of interchangeable phantom inserts with different pathologies
- 3 | Realistic view of female external genital
- 4 | Excellent ultrasound image quality

- | Handling and manipulation of transvaginal and transabdominal transducers
- | Interpretation of sonogram
- Visualization and localization of anatomies and pathologies

Three variations of training sets. Select the one which best meet your requirements. ANATOMY



DESCRIPTIONS			
SPECIFICATIONS		MATERIALS	REPLACEMENT PARTS
Torso size: W34 x D33 x H24cm W13.5 x D13 x H9.5 inch SET INCLUDES	weight: 7.2kg / 15.9 lbs	Soft resin, Polyurethane elastomer Latex free	41909-030 ultrasound pathology unit 41909-040 ultrasound ectopic pregnancy unit 41909-050 ultrasound first trimester pregnancy unit
US-10 Female Pelvic Ultrasound Phantom	US-10A Female Pelvic Ultrasound Phantom	US-10B Obstetrics Pelvic Ultrasound Phantom	41909-060 skin for US-10 / 10A / 10B
1 lower torso manikin	1 lower torso manikin	1 lower torso manikin	
1 ultrasound pathology unit	1 ultrasound pathology unit	1 ultrasound first trimester pregnancy ur	nit
1 ultrasound ectopic pregnancy unit	1 ultrasound first trimester pregnancy unit	1 ultrasound ectopic pregnancy unit	
1 talcum powder	1 ultrasound ectopic pregnancy unit	1 talcum powder	
1 Syringe	1 talcum powder	1 syringe	
1 storage case	1 syringe	1 storage case	

- 1 instruction manual
- KYOTO KAGAKU 🔿 KYOTOKAGAKU AMERICA INC.

2 storage cases

instruction manual

6

1 instruction manual

https://www.kyotokagaku.com/en/ Q

* Ultrasound scanner is not included.

