

# Adult CRiSis™ Manikin®

Product Name	<b>Adult CRiSis™ Manikin®</b>
Addresses Training Needs/ Curricula	<ul style="list-style-type: none"> <li>• Realistic training and practice of airway management and CPR education</li> </ul>
Configuration	<ul style="list-style-type: none"> <li>• Simulating a full-body adult. It combines all skill stations on 1 manikin.</li> <li>• Realistic anatomy and landmarks.</li> <li>• Upgradeable to include interactive ECG.</li> <li>• Advanced Airway Management, Auscultation, CPR monitoring and intraosseous infusion.</li> </ul>
Dimension	<ul style="list-style-type: none"> <li>• Length: 5'5" (65 inches)</li> <li>• Ship weight: 88 lbs</li> </ul>
Product Code	LF03890 Adult CRiSis™ Manikin® LF01144 Auscultation Smartscope™ LF01148 Auscultation Remote LF03406 CPR Metrix with iPad LF03670 ECG Simulator LF03684 Advanced Airway Management Head LF03840 Articulating Legs with Intraosseous Infusion and IM Injection



# Be READY



**Nasco Healthcare Inc.**  
 16 Simulaids Drive, Saugerties, NY, 12477 USA  
 Phone: 1-833-NASCOHC (627-2642)  
 Email: info@nascohealthcare.com  
 www.nascohealthcare.com



Ready to make a difference.  
Ready to save lives.

## Adult CRiSis™ Manikin®

The most comprehensive, durable, low-cost EMS training simulator on the market

- All-in-one complete ACLS/ACS training capabilities in multiple configurations; full body or torso
- Light weight
- Low maintenance
- Ready-to-use
- Durable construction and airway
- 5 years warranty



Complete CRiSis™ Manikin® (LF03890) Sh. wt. 88 lbs.

Deluxe CRiSis™ Manikin® Torso (LF03958U) includes an ECG arrhythmia simulator.

Deluxe CRiSis™ Manikin® Torso with Advanced Airway Management (LF03983U) includes laryngospasms (manual) and tongue edema (manual).

Advanced "Airway Larry" Torso with Defibrillation Features (LF03960U) includes tongue edema and laryngospasm.

## Complete CRiSis™ Manikin®

Realistic, comprehensive emergency training

Enables training on airway management and proper CPR techniques, including:

**Constricted and Dilated Eyes**

**CPR:** Jaw Thrust; palpable and visual landmarks; fully articulated head, neck, and jaw; compressible chest for adult and child.

**Bilateral, Manual Carotid Pulse**

**Airway Management:** Oral, nasal, and digital intubation; suctioning capabilities; practice of Sellick maneuver; realistic anatomy. Compatible airway devices include E.T., E.O.A., E.T.L, L.M.A, E.G.T.A., Combitube® and King systems.

**Blood Pressure Arm:** Auscultatory BP with 5 Korotkoff sounds and palpable radial pulse that can be turned on and off.

**Auscultation:** 6 anterior heart sites with 12 conditions; 5 anterior, 2 mid axillary lung sites with 16 conditions.

**Defibrillation Chest:** Manual, semi-automatic, automatic defibrillation; compatible with all standard brands and types of defibrillators; internal load box absorbs upto 360 joules of every shock.

**4-Lead ECG**

**I.V. Arm:** Articulated at the biceps; intramuscular injections; realistic flash-back.

This manikin is intended to be used as an instructional aid in the teaching and demonstrating of live rescue training. Actual product may vary slightly from photo. Nasco Healthcare reserves the right to change product color, materials, supplies, or function as needed.

